

Miracle Formula for Cleaning Wood Furniture

3 parts Boiled Linseed Oil
2 parts Turpentine
1 part vinegar

For a normal batch, I use 'third cups' for measuring (1 cup Linseed, 2/3 cup turpentine, 1/3 cup vinegar)

Using a wire whisk, thoroughly mix the three items until the mixture resembles Buttermilk. This may take up to 5 minutes of briskly whisking! Using a bristle brush, liberally apply the product to a manageable working area, rubbing the product in a crevices and surfaces. Let stand for about 5 minutes, working the product in again and again with your brush. Remove residue with a cotton cloth.

Repeat until product removed is clean and surface is restored.

After your piece is cleaned, use your wipe-off rag and lightly go over the piece, leaving behind a VERY thin layer of product. Let dry for at least a week.

Discard any un-used material after 5 days.

© 2002-2004 Southeast Antique Radio Society. ALL RIGHTS RESERVED. UNAUTHORIZED DUPLICATION PROHIBITED.

Caution: Performing repairs on radios could be dangerous. SARS assumes no responsibility for accidents resulting from any information contained in this web site.